



## DEVICE DETAILS

NAME OF DEVICE	Disposable Medical 3ply Face Mask
ESTABLISHMENT NAME	NAEL MARKETING
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA6541621-75062
BRAND NAME	QJMDM ( QIANJIANG KINGPHAR MEDICAL MATERIAL)
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	A 3ply (three layers) medical grade face mask made from a melt-blown material placed between non-woven fabric. The melt-blown material acts as the filter that stops microbes from entering or exiting the mask. Most medical grade masks feature pleats or folds.
DEVICE INTENDED PURPOSE	Medical Face Mask with efficient protection to protect user from air-borne germs and viruses. It also prevents viruses sneezed from user's nose to spread outside
VALIDITY DATE OF REGISTRATION	08/10/2021 - 07/10/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Disposable Medical 3ply Face Mask	QJMDM TYPE II R