



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	4-Ply Premium Face Mask
<b>ESTABLISHMENT NAME</b>	HIL MEDIC SDN. BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA5329521-76076
<b>BRAND NAME</b>	HIL Medic
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	It is a Non-woven Mask which meets all current ASTM F2100 level 3 and EN 14683 IIR standards for face masks. Used to cover the mouth and nose. They serve as a barrier for contaminants from getting released from the mouth of a sick person and spread to the environment, it suitable for use for surgical procedures purpose or by people in the public
<b>DEVICE INTENDED PURPOSE</b>	To be used to protect both the users from transfer of microorganisms, body fluids, and particulate material.
<b>VALIDITY DATE OF REGISTRATION</b>	15/10/2021 - 14/10/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	4 Ply Premium Face Mask	FM005/01 - Colours
2	4 Ply Premium Face Mask	FM005/02 - Pattern