



DEVICE DETAILS

NAME OF DEVICE	3 Ply Medical Face Mask
ESTABLISHMENT NAME	MEDTECH & HEALTH SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA10408421-76570
BRAND NAME	Medtech & Health Sdn Bhd
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	3 Ply Medical Face Mask protects against moisture from the outside and protects against 98% of dust, pollen, bacteria, viruses and other airborne particles. Outer layer = 25g/m2 PP Brown Spunbond Filter layer = 25g/m2 meltblown, 98% BFE Inner layer = ES thermobonded soft Plastic and wire nosepiece 3mm latex-free round elastic earloop
DEVICE INTENDED PURPOSE	a flexible, loose-fitting mask designed to be placed over the mouth and/or nose and chin fitted with the head harness which can be head or ears attachment of a wearer to permit normal breathing while protecting the wearer from the transfer of particles from the environment.
VALIDITY DATE OF REGISTRATION	21/10/2021 - 20/10/2026

LIST OF DEVICE



NO	NAME OF DEVICE	IDENTIFIER
1	3 Ply Medical Face Mask	09551009700039
2	3 Ply Medical Face Mask	09551009700022
3	3 Ply Medical Face Mask	09551009700015
4	3 Ply Medical Face Mask	09551009700008