



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	3 Ply Medical Face Mask
<b>ESTABLISHMENT NAME</b>	MEDTECH & HEALTH SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA10408421-76570
<b>BRAND NAME</b>	Medtech & Health Sdn Bhd
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	3 Ply Medical Face Mask protects against moisture from the outside and protects against 98% of dust, pollen, bacteria, viruses and other airborne particles. Outer layer = 25g/m2 PP Brown Spunbond Filter layer = 25g/m2 meltblown, 98% BFE Inner layer = ES thermobonded soft Plastic and wire nosepiece 3mm latex-free round elastic earloop
<b>DEVICE INTENDED PURPOSE</b>	a flexible, loose-fitting mask designed to be placed over the mouth and/or nose and chin fitted with the head harness which can be head or ears attachment of a wearer to permit normal breathing while protecting the wearer from the transfer of particles from the environment.
<b>VALIDITY DATE OF REGISTRATION</b>	21/10/2021 - 20/10/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	3 Ply Medical Face Mask	09551009700039
2	3 Ply Medical Face Mask	09551009700022
3	3 Ply Medical Face Mask	09551009700015
4	3 Ply Medical Face Mask	09551009700008