



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	SURGICAL FACE MASK
<b>ESTABLISHMENT NAME</b>	GLOBAL GIFUTO SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GA9665921-77096
<b>BRAND NAME</b>	BL/WA-FM
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	Surgical mask is composed of mask body, nose clip, and mask belt. The mask body is composed of spunbonded non-woven fabric and meltblown filter material. The nose clip is made of polyvinyl chloride (PVC) wrapped metal wire
<b>DEVICE INTENDED PURPOSE</b>	it is used to cover the mouth, nose and jaw of the user, and to provide a physical barrier to prevent the direct penetration of pathogenic microorganisms, body fluids, and particles
<b>VALIDITY DATE OF REGISTRATION</b>	01/11/2021 - 31/10/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	SURGICAL FACE MASK	WA-FM