



DEVICE DETAILS

NAME OF DEVICE	SURGICAL FACE MASK
ESTABLISHMENT NAME	GLOBAL GIFUTO SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA9665921-77096
BRAND NAME	BL/WA-FM
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	Surgical mask is composed of mask body, nose clip, and mask belt. The mask body is composed of spunbonded non-woven fabric and meltblown filter material. The nose clip is made of polyvinyl chloride (PVC) wrapped metal wire
DEVICE INTENDED PURPOSE	it is used to cover the mouth, nose and jaw of the user, and to provide a physical barrier to prevent the direct penetration of pathogenic microorganisms, body fluids, and particles
VALIDITY DATE OF REGISTRATION	01/11/2021 - 31/10/2026

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	SURGICAL FACE MASK	WA-FM