



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	MEDICAL FACE MASK 4 LAYERS
<b>ESTABLISHMENT NAME</b>	SIMPLY K SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA4606921-78010
<b>BRAND NAME</b>	SIMPLY K
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	4-PLY PREMIUM MEDICAL FACE MASK WITH DOUBLE MELT-BLOWN PROTECTION. EASY TO BREATHE THROUGH. ODORLESS AND NON-IRRITATING. HIGH FILTRATION CAPACITY. HELPS FILTER OUT DUST, BACTERIA, SMOKE & POLLEN.
<b>DEVICE INTENDED PURPOSE</b>	Medical face mask that covers the user's nose and mouth and provides a physical barrier from transmission of viruses, microorganisms and dust. It is intended to be use by healthcare provider and other persons for medical purposes and to control infectious diseases, particularly in epidemic or pandemic situations.
<b>VALIDITY DATE OF REGISTRATION</b>	17/11/2021 - 16/11/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	MEDICAL FACE MASK 4 LAYERS	MEDICAL FACE MASK 4 LAYERS ADULT - EAR LOOP
2	MEDICAL FACE MASK 4 LAYERS	MEDICAL FACE MASK 4 LAYERS ADULT- HEAD LOOP
3	MEDICAL FACE MASK 4 LAYERS	MEDICAL FACE MASK 4 LAYERS CHILDREN - EAR LOOP
4	MEDICAL FACE MASK 4 LAYERS	MEDICAL FACE MASK 4 LAYERS CHILDREN - HEAD LOOP