



DEVICE DETAILS

NAME OF DEVICE	SURGICAL FACE MASK
ESTABLISHMENT NAME	ICON PACKAGING SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA2221021-78544
BRAND NAME	SUNWAY
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	A single use non sterile and latex free device made of 2 layer of non woven and 1 layer of melt blown. Adjustable to fit faces of various shapes and sizes. The middle layer of face mask (Melt blown) which tested for Bacterial Filtration Efficiency that has sufficient protection against haze particles, dust and harmful droplets. It also high breathability, hypoallergenic, soft and lightweight. This device is appropriate to be use by healthcare professional during surgical procedures to protect both the patient and healthcare personnel.
DEVICE INTENDED PURPOSE	A Surgical face mask is meant to help block large particle droplets, splashes or splatter that may contain germs (viruses & bacteria), keeping it from reaching to your mouth.
VALIDITY DATE OF REGISTRATION	19/11/2021 - 18/11/2026

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	3 PLY SURGICAL FACE MASK	SW001
2	3 PLY SURGICAL FACE MASK	SW003
3	3 PLY SURGICAL FACE MASK	SW006