



DEVICE DETAILS

NAME OF DEVICE	HANDHELD HANDHELD PULSE OXIMETER
ESTABLISHMENT NAME	IHB HEALTHCARE SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GC5337721-78806
BRAND NAME	AEON TECHNOLOGY
MEDICAL DEVICE CATEGORY	MD 1302 - Monitoring devices of vital physiological parameters
DEVICE GROUPING TYPE	SYSTEM
DEVICE DESCRIPTION	<p>Aeon Technology A36 Handheld Handheld Pulse Oximeter is a is a non-invasive, handheld patient Oximeter. It operates on alkaline or rechargeable battery power supply. It is compact, small, light, and easy for learning and handling. It is suitable for monitoring adult and child patients. It is widely used only in Professional medical organizations, such as, hospitals, clinics.</p> <p>Parameters measured by the oximeter include: arterial oxygen saturation (SpO2), pulse rate (PR), bargraph and plethysmogram. The oximeter measures these parameters through a SpO2 sensor and displays them on the color TFT LCD screen after certain further processing. The oximeter is operated and controlled by the buttons on the front panel. It adopts a 2.8 inch color TFT LCD screen in displaying measurements and in supplementary status indication.</p>
DEVICE INTENDED PURPOSE	Aeon Technology Handheld Pulse Oximeter is a non-invasive device intended for monitoring adult and child patients. It is widely used only in Professional medical organizations, such as, hospitals, clinics.
VALIDITY DATE OF REGISTRATION	21/11/2021 - 20/11/2026

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Handheld Pulse Oximeter	A36
2	Adult finger SpO2 sensor	N.A.
3	USB line	N.A.
4	1.5V Alkaline Battery	N.A.
5	Optical disk	N.A.
6	Neonatal finger SpO2 sensor	N.A.
7	Child finger SpO2 sensor	N.A.

NO	NAME OF DEVICE	IDENTIFIER
8	Temperature probe	N.A.
9	Rechargeable Li battery	N.A.



MEDICAL DEVICE AUTHORITY