



## DEVICE DETAILS

|                                      |  |
|--------------------------------------|--|
| <b>NAME OF DEVICE</b>                | KN95 PREMIUM PROTECTIVE SURGICAL FACE MASK   |
| <b>ESTABLISHMENT NAME</b>            | WORLDTREND ENTERPRISE  |
| <b>ROLE OF ESTABLISHMENT</b>         | AUTHORIZED REPRESENTATIVE  |
| <b>REGISTRATION NO</b>               | GA6535221-79322  |
| <b>BRAND NAME</b>                    | MEDI-CARE  |
| <b>MEDICAL DEVICE CATEGORY</b>       | SINGLE-USE DEVICES   |
| <b>DEVICE GROUPING TYPE</b>          | SINGLE   |
| <b>DEVICE DESCRIPTION</b>            | Disposable surgical face masks embedded with nose clip and ear loops or with ties. This surgical face mask made of 5-ply of premium non-woven fabric with Bacterial Filtration Efficiency (BFE) >99%. Its involved in 4-fold design for comfortable fit with strong and elastic ear loops.   |
| <b>DEVICE INTENDED PURPOSE</b>       | This product is used by medical staff for disposable health care in normal environment. Disposable surgical face masks embedded with nose clip and ear loops or with ties. This surgical face mask made of 5-ply of premium non-woven fabric with Bacterial Filtration Efficiency (BFE) >99%. Its involved in 4-fold design for comfortable fit with strong and elastic ear loops. |
| <b>VALIDITY DATE OF REGISTRATION</b> | 26/11/2021 - 25/11/2026  |

## LIST OF DEVICE

| NO | NAME OF DEVICE                             | IDENTIFIER |
|----|--|------------|
| 1  | KN95 PREMIUM PROTECTIVE SURGICAL FACE MASK | MC-004     |