



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	FINGERTIP PULSE OXIMETER
<b>ESTABLISHMENT NAME</b>	TST HEALTH CARE SDN. BHD.
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GC3239621-80138
<b>BRAND NAME</b>	SONOSAT F02W
<b>MEDICAL DEVICE CATEGORY</b>	MD 1302 - Monitoring devices of vital physiological parameters
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	FINGERTIP PULSE OXIMETER IS A PORTABLE NON-INVASIVE SPOT-CHECK MEASUREMENT DEVICE APPLICABLE TO CHECK BLOOD OXYGEN SATURATION (SPO <sub>2</sub> ), PULSE MONITORING (BPM), AND BLOOD PERFUSION INDEX (PI) THROUGH A PATIENT'S FINGER. THIS MEDICAL DEVICE CAN BE REUSED. NOT FOR CONTINUING MONITORING.
<b>DEVICE INTENDED PURPOSE</b>	TO BE USED FOR HOME CARE AND MEDICAL OUTPATIENT DEPARTMENT MEASURE PULSE OXYGEN SATURATION AND PULSE RATE OF ADULTS AND PAEDIATRICS.
<b>VALIDITY DATE OF REGISTRATION</b>	01/12/2021 - 30/11/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	FINGERTIP PULSE OXIMETER	SONOSAT F02W