



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	KF94 PREMIUM PROTECTIVE SURGICAL FACE MASK
<b>ESTABLISHMENT NAME</b>	WORLDTREND ENTERPRISE
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GA3196221-80190
<b>BRAND NAME</b>	MEDI-CARE
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	Disposable surgical face masks embedded with nose clip and ear loops or with ties. This surgical face mask made of 4-ply of premium non-woven fabric with Bacterial Filtration Efficiency (BFE) >99% and Particle Filtration Efficiency (PFE) >95%. Its involved in 4-fold design for comfortable fit with strong and elastic ear loops.
<b>DEVICE INTENDED PURPOSE</b>	Intended for single use by operating room personnel and other general healthcare workers to protect both patients and healthcare workers against transfer of microorganisms, blood and body fluids, and airborne particulates.
<b>VALIDITY DATE OF REGISTRATION</b>	05/12/2021 - 04/12/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	KF94 PREMIUM PROTECTIVE SURGICAL FACE MASK	MC-003