

## DEVICE DETAILS

NAME OF DEVICE	KF94 PREMIUM PROTECTIVE SURGICAL FACE MASK
ESTABLISHMENT NAME	WORLDTREND ENTERPRISE
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA3196221-80190
BRAND NAME	MEDI-CARE
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	Disposable surgical face masks embedded with nose clip and ear loops or with ties. This surgical face mask made of 4-ply of premium non-woven fabric with Bacterial Filtration Efficiency (BFE) >99% and Particle Filtration Efficiency (PFE) >95%. Its involved in 4-fold design for comfortable fit with strong and elastic ear loops.
DEVICE INTENDED PURPOSE	Intended for single use by operating room personnel and other general healthcare workers to protect both patients and healthcare workers against transfer of microorganisms, blood and body fluids, and airborne particulates.
VALIDITY DATE OF REGISTRATION	05/12/2021 - 04/12/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	KF94 PREMIUM PROTECTIVE SURGICAL FACE MASK	MC-003