



DEVICE DETAILS

NAME OF DEVICE	Antimicrobial KF94 Mask
ESTABLISHMENT NAME	BRAVOMED SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA7881621-80490
BRAND NAME	EG GUARD
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Antimicrobial KF94 Mask are made of 4-ply of water resistant outer layer, non woven filter to filter out the dust and germs, excellent filtering melt blown fabric with electrostatic charge, soft and fine inner non woven skin friendly. They have bacterial filtration efficiency and particles filtration efficiency for use by humans in any situation such as pandemic, medical sectors and common use for health concern.
DEVICE INTENDED PURPOSE	Antimicrobial KF94 Mask is intended to be worn by Medical Personnel during surgical or other Medical Procedures to protect both the patient and the operating personnel from transfer of microorganisms, body fluid and particulate material transfer.
VALIDITY DATE OF REGISTRATION	06/12/2021 - 05/12/2026

LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

NO	NAME OF DEVICE	IDENTIFIER
1	Antimicrobial 2D Mask	KF94 2D MASK
2	Antimicrobial 3D Mask	KF94 3D MASK