



DEVICE DETAILS

| | |
|--------------------------------------|--|
| NAME OF DEVICE | FINGERTIP PULSE OXIMETER |
| ESTABLISHMENT NAME | MALAYSIAN HEALTHCARE SDN BHD |
| ROLE OF ESTABLISHMENT | AUTHORIZED REPRESENTATIVE |
| REGISTRATION NO | GC2345922-85272 |
| BRAND NAME | YUWELL |
| MEDICAL DEVICE CATEGORY | MD 1302 - Monitoring devices of vital physiological parameters |
| DEVICE GROUPING TYPE | FAMILY |
| DEVICE DESCRIPTION | The finger pulse oximeter features in small volume, low power consumption, convenient operation and portable. It is only necessary for patient to put one of his fingers into a fingertip photoelectric sensor for measurement, and then the screen will display the measured value of oxyhemoglobin saturation. It has been proved in clinical experiments that it features in rather than high precise and repeatability |
| DEVICE INTENDED PURPOSE | For hospital and homecare to supply oxygen to patients. Supplying oxygen to patients can assist in the treatment of heart head blood-vessel, respiratory system, such as chronic obstructive pneumonia disease, rehabilitation and hypoxia conditions |
| VALIDITY DATE OF REGISTRATION | 24/01/2022 - 23/01/2027 |

LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

| NO | NAME OF DEVICE | IDENTIFIER |
|----|--------------------------|------------|
| 1 | Fingertip Pulse Oximeter | YX306 |
| 2 | Fingertip Pulse Oximeter | YX102 |