



DEVICE DETAILS

NAME OF DEVICE	LEXCO 3D MEDICAL FACE MASK
ESTABLISHMENT NAME	STARWIL (M) SDN. BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA9207322-86431
BRAND NAME	LEXCO
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Face masks are made of 3-ply of water resistant outer layer, non woven filter to filter out the dust and germs, excellent filtering melt blown fabric with electrostatic charge, soft and fine inner non woven skin friendly and suitable for sensitive skin. They have bacterial filtration efficiency and particles filtration efficiency for use by humans in any situation such as pandemic, medical sectors and common use for health concern. A new 3D design for easy talking while using the face mask.
DEVICE INTENDED PURPOSE	3 layers filtering face mask respirator and used as to reduce or prevent the transmission of disease. It's single use product and is not recommended for multiple use.
VALIDITY DATE OF REGISTRATION	25/02/2022 - 24/02/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	LEXCO 3D MEDICAL FACE MASK	3D-B
2	LEXCO 3D MEDICAL FACE MASK	3D-W