



## DEVICE DETAILS

NAME OF DEVICE	KF94 4PLY MEDICAL FACE MASK
ESTABLISHMENT NAME	CONCEPT HERO SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA1224222-88199
BRAND NAME	AMDPRO+
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	The 4 ply surgical mask with an added extra layer with one more filtering layer. They also have adjustable nose strip to give maximum protection and comfort to the user. With an added extra layer ensuring that you remain protected with 98% bacterial and particle filtration. A melt blown layer is critical as a filter because it provides a layer that contains droplets by electrostatically absorbing them on the surface. This makes it so that dust, particles, and bacteria cannot infiltrate the mask.
DEVICE INTENDED PURPOSE	Medical/Surgical mask is meant to help block large-particle droplet, splashes, sprays or splatter that contain germs (viruses and bacteria), keeping it from reaching your mouth and nose. Surgical mask may also help reduce exposure of your saliva n respiratory secretions to others
VALIDITY DATE OF REGISTRATION	18/03/2022 - 17/03/2027

## LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

NO	NAME OF DEVICE	IDENTIFIER
1	KF94 DISPOSABLE MEDICAL FACE MASK, 4PLY	FM-AP+KF94-K-M-200
2	KF94 DISPOSABLE MEDICAL FACE MASK, 4PLY	FM-AP+KF94-M-220