



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	MOVEN Lightweight Transport Wheelchair
<b>ESTABLISHMENT NAME</b>	BEST MADE PHARMACEUTICAL SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GA9618422-91832
<b>BRAND NAME</b>	Lightweight Transport Wheelchair / MOVEN
<b>MEDICAL DEVICE CATEGORY</b>	ASSISTIVE PRODUCTS FOR PERSONS WITH DISABILITY
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	Lightweight wheelchair made of aluminum frame, foldable backrest, fixed armrest, detachable legrest, double crossbar and 8"/12" front/back wheels.
<b>DEVICE INTENDED PURPOSE</b>	To help individuals who have impairments that limit their ability to walk to move.
<b>VALIDITY DATE OF REGISTRATION</b>	07/05/2022 - 06/05/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	MOVEN Lightweight Transport Wheelchair	BM01LABJ-12