




DEVICE DETAILS

NAME OF DEVICE	MOVEN Lightweight Transport Wheelchair
ESTABLISHMENT NAME	BEST MADE PHARMACEUTICAL SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA9618422-91832
BRAND NAME	Lightweight Transport Wheelchair / MOVEN
MEDICAL DEVICE CATEGORY	ASSISTIVE PRODUCTS FOR PERSONS WITH DISABILITY
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	Lightweight wheelchair made of aluminum frame, foldable backrest, fixed armrest, detachable legrest, double crossbar and 8"/12" front/back wheels.
DEVICE INTENDED PURPOSE	To help individuals who have impairments that limit their ability to walk to move.
VALIDITY DATE OF REGISTRATION	07/05/2022 - 06/05/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	MOVEN Lightweight Transport Wheelchair 	BM01LABJ-12