



DEVICE DETAILS

NAME OF DEVICE	Disposable Medical Face Mask
ESTABLISHMENT NAME	ALLIANCE INTEGRITY TRADING SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA4172222-82952
BRAND NAME	MILKON
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Disposable medical face mask this product is composed of non-woven cloth, filter cloth, nose clip, mask belt (rubber band); The mask is made of two layers of non-woven cloth and one layer of filter cloth by folding ultrasonic welding. It is non-sterile and disposable. The filtration efficiency of BFE bacteria was $\geq 98\%$, the total number of bacterial colonies was ≤ 100 cfu/g, and the PH value was 6.0-8.0.
DEVICE INTENDED PURPOSE	This mask is a disposable product, suitable for the health care of the wearer in the general medical environment without the risk of body fluids and spillage and general health care in public health places.
VALIDITY DATE OF REGISTRATION	14/01/2022 - 13/01/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Disposable Medical Face Mask	MLK-05
2	Disposable Medical Face Mask	MLK-05