



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	FINGERTIP OXIMETER
<b>ESTABLISHMENT NAME</b>	UCREST BERHAD
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GC8852522-92820
<b>BRAND NAME</b>	CREATIVE
<b>MEDICAL DEVICE CATEGORY</b>	MD 1302 - Monitoring devices of vital physiological parameters
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	Fingertip Oximeter (PC-60FW) is a non-invasive, non-sterile portable device intended for measuring the pulse rate and functional oxygen saturation (SpO <sub>2</sub> ) through a patient's finger. It is applicable for spot checking of SpO <sub>2</sub> level and pulse rate of adult and paediatric patients at homes and medical clinics settings. It is lightweight, small in size, low power consumption, convenient operation and portability. It only required patient to place one finger into fingertip photoelectric sensor for diagnosis, and the display screen will directly show measured value of functional oxygen saturation and pulse rate.
<b>DEVICE INTENDED PURPOSE</b>	Fingertip Oximeter is a non-invasive, non-sterile portable device intended for measuring the pulse rate and functional oxygen saturation (SpO <sub>2</sub> ) through a patient's finger. It is indicated for spot checking of SpO <sub>2</sub> level and pulse rate of adult and paediatric patients at homes and medical clinics settings.
<b>VALIDITY DATE OF REGISTRATION</b>	13/05/2022 - 12/05/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	FINGERTIP OXIMETER	11000207