



DEVICE DETAILS

NAME OF DEVICE	MACare 3-ply Medical Disposable Face Mask
ESTABLISHMENT NAME	MATC TECHNOLOGY (M) SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA4855522-93990
BRAND NAME	MACare
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Face mask outer layer is made with non-woven fabric; mid-layer is made with non-woven melt brown fabric which provides >95% filtration efficiency, inner layer is made with sky-friendly non-woven fabric. Ear-lop is made with latex free material. And nose strip is from PP single core 3mm. Single time use face mask
DEVICE INTENDED PURPOSE	The mask is intended to be worn by medical personnel during surgical or other medical procedures to protect both the patient and the operating personnel from transfer of microorganisms, body fluid and particulate material
VALIDITY DATE OF REGISTRATION	08/06/2022 - 07/06/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	3-ply Disposable Medical Facemask (Midnight Sky)	553931-69700080
2	3-ply Disposable Medical Facemask (Lemon Rock)	553931-69700090
3	3-ply Disposable Medical Facemask (Fairy Floss)	553931-69700100
4	3-ply Disposable Medical Facemask (Mistry Leaf)	553931-69700110
5	3-ply Disposable Medical Facemask (Couple Pegar)	553931-69700120
6	3-ply Disposable Medical Facemask (Bunga Raya)	553931-69700130