



DEVICE DETAILS

NAME OF DEVICE	N95 SMP ANTIBACTERIAL MEDICAL FACE MASK
ESTABLISHMENT NAME	DEXINA FLORAPRO RESOURCES (M) SDN. BHD.
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA3139322-94046
BRAND NAME	TRUVIGOR
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	TRADE MARK : TRUVIGOR MODEL : N95 SMP ANTIBACTERIAL MEDICAL FACE MASK RATING EN149 CLASSIFICATION =FFP2 SPLASH RESISTANCE TYPE :IIR [FACE MASK TYPE RESPIRATORY PROTECTIVE DEVICE TO BE USED IN AND/OR OPERATING THEATRE OTHER MEDICAL SETTINGS] TEST REPORT :2021CE1949 , 2022CE0104, 03-14B-02-2022-39-03-1 , 2022CE010-S1 , 03-14B-02-2022-39-04-2 PRODUCT STANDARD: EN14683:2019 PRODUCT : MEDICAL FACE MASK SHAPE: 3D PROTECTIVE 5 PLY FACE MASK AVAILABILITY :SINGLE USE ONLY COUNTRY OF ORIGIN :MALAYSIA (MADE IN MALAYSIA) MANUFACTURER: DEXINA FLORAPRO RESOURCES (M) SDN BHD FEATURE :ANTIBACTERIAL ,ANTIVIRAL ,DEODORIZATION ,SAFE,COMFORTABLE ,STYLISH
DEVICE INTENDED PURPOSE	TO BE WORN BY USER AND ACT AS PROTECTIVE BARRIER FOR NOSE AND MOUTH FROM LARGE PARTICLE DROPLETS SPLASHES , SPRAYS OR SPLATTER THAT MAY CONTAIN VIRUSES AND BACTERIA OR ANY OTHER AIR - BORNE PARTICULATES .THIS FACE MASK USED AS A SURGICAL MASK .
VALIDITY DATE OF REGISTRATION	08/06/2022 - 07/06/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	N95 SMP ANTIBACTERIAL MEDICAL FACE MASK	N95 SMP