



DEVICE DETAILS

NAME OF DEVICE	FACE MASK MEDICAL DISPOSABLE
ESTABLISHMENT NAME	MEDIRIAL SDN. BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA7973622-94186
BRAND NAME	M-CURE
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	4ply face mask intended to be use at hospital, dental, clean room, food preparation and industrial environments. no-irritating to skin and permits air passage. comes with 4 layer. Out layer, Blue 20 polypropylene spun bonded, middle layer, blue 20g polypropylene spun bonded, middle layer, white 25g polypropylene spun bonded, inner layer, white 25g polypropylene
DEVICE INTENDED PURPOSE	4ply face mask intended to be use at hospital, dental, clean room, food preparation and industrial environments. no-irritating to skin and permits air passage.
VALIDITY DATE OF REGISTRATION	08/06/2022 - 07/06/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	FACE MASK MEDICAL DISPOSABLE	087400-(W)
2	FACE MASK MEDICAL DISPOSABLE	087400-(B)
3	FACE MASK MEDICAL DISPOSABLE	087400-(G)
4	FACE MASK MEDICAL DISPOSABLE	087400-(PR)
5	FACE MASK MEDICAL DISPOSABLE	087400-(PK)
6	FACE MASK MEDICAL DISPOSABLE	088400-(W)
7	FACE MASK MEDICAL DISPOSABLE	088400-(B)
8	FACE MASK MEDICAL DISPOSABLE	088400-(G)
9	FACE MASK MEDICAL DISPOSABLE	088400-(PR)
10	FACE MASK MEDICAL DISPOSABLE	088400-(PK)

NO	NAME OF DEVICE	IDENTIFIER
11	FACE MASK MEDICAL DISPOSABLE	089400-(W)
12	FACE MASK MEDICAL DISPOSABLE	089400-(B)
13	FACE MASK MEDICAL DISPOSABLE	089400-(G)
14	FACE MASK MEDICAL DISPOSABLE	089400-(PR)
15	FACE MASK MEDICAL DISPOSABLE	089400-(PK)



MEDICAL DEVICE AUTHORITY