



DEVICE DETAILS

NAME OF DEVICE	3-Ply Disposable Single-Use Surgical Face Mask
ESTABLISHMENT NAME	JAIHIND ENGINEERING SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA6750222-94872
BRAND NAME	JH 3-Ply Disposable Single-Use Surgical Face Mask
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	A disposable single-use surgical face mask made up of 3 layers. First/outer layer: polypropylene spunbond nonwoven fabric, middle/filter layer: polypropylene melt-blown fabric that has BFE > 99% and third/inner layer: polypropylene spunbond nonwoven fabric.
DEVICE INTENDED PURPOSE	To protect user against or to reduce the risks of the transfer of microorganisms, blood and bodily fluids and airborne particles.
VALIDITY DATE OF REGISTRATION	17/06/2022 - 16/06/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	3-Ply Disposable Single-Use Surgical Face Mask	EARLOOP PREMIUM PROCEDURE 3-PLY FACE MASK