



DEVICE DETAILS

| | |
|-------------------------------|---|
| NAME OF DEVICE | 3-Ply Disposable Single-Use Surgical Face Mask |
| ESTABLISHMENT NAME | JAIHIND ENGINEERING SDN BHD |
| ROLE OF ESTABLISHMENT | MANUFACTURER |
| REGISTRATION NO | GA6750222-94872 |
| BRAND NAME | JH 3-Ply Disposable Single-Use Surgical Face Mask |
| MEDICAL DEVICE CATEGORY | SINGLE-USE DEVICES |
| DEVICE GROUPING TYPE | SINGLE |
| DEVICE DESCRIPTION | A disposable single-use surgical face mask made up of 3 layers. First/outer layer: polypropylene spunbond nonwoven fabric, middle/filter layer: polypropylene melt-blown fabric that has BFE > 99% and third/inner layer: polypropylene spunbond nonwoven fabric. |
| DEVICE INTENDED PURPOSE | To protect user against or to reduce the risks of the transfer of microorganisms, blood and bodily fluids and airborne particles. |
| VALIDITY DATE OF REGISTRATION | 17/06/2022 - 16/06/2027 |

LIST OF DEVICE

| NO | NAME OF DEVICE | IDENTIFIER |
|----|--|---|
| 1 | 3-Ply Disposable Single-Use Surgical Face Mask | EARLOOP PREMIUM PROCEDURE 3-PLY FACE MASK |