



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	Surgical Face Mask
<b>ESTABLISHMENT NAME</b>	U CARE PROTECTION INDUSTRIES SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA3161022-94988
<b>BRAND NAME</b>	Surgical Face Mask / Daxxo
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	A single use non sterile and latex free device made of 2 layer of non woven and 1 layer of melt blown. Adjustable to fit faces of various shapes and sizes. The middle layer of face mask (Melt blown) which tested for Bacterial Filtration Efficiency that has sufficient protection against haze particles, dust and harmful droplets. It also high breathability, hypoallergenic, soft and lightweight. This device is appropriate to be use by healthcare professional during surgical procedures to protect both the patient and healthcare personnel.
<b>DEVICE INTENDED PURPOSE</b>	A disposable single use surgical face mask is to be used as a barrier to prevent air-borne transfer of infectious, large particle droplets and various micro-organisms. It is designed to be used in hospital, clinic, operating theatre, dental, cleanroom and critical environment.
<b>VALIDITY DATE OF REGISTRATION</b>	20/06/2022 - 19/06/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Surgical Face Mask	DX4P101