



## DEVICE DETAILS

NAME OF DEVICE	3ply Surgical Face Mask
ESTABLISHMENT NAME	IMEDIC SUPPLY SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA3623222-98283
BRAND NAME	I?Mask
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Three-ply construction for excellent particle and bacterial filtration efficiency. These masks are available with spandex ear loops or polypropylene ties, ultrasonically welded to maintain softness and protect against particulates. An encased, 12cm (4¾") malleable steel noseband creates a secure facial seal and improves wearer comfort. Latex-free. Packaged 50 masks per pack, 10 packs per case. Mask dimensions: 17.5L x 9.5Hcm
DEVICE INTENDED PURPOSE	The Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non-sterile.
VALIDITY DATE OF REGISTRATION	26/07/2022 - 25/07/2027



## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	3 PLY SURGICAL FACEMASK (EARLOOP) COLOR: BLUE	IMBEBFE99
2	3 PLY SURGICAL FACEMASK (EARLOOP) COLOR: WHITE	IMWEBFE99
3	3 PLY SURGICAL FACEMASK (HEADLOOP) COLOR: BLUE	IMHBBFE99
4	3 PLY SURGICAL FACEMASK (HEADLOOP) COLOR: WHITE	IMHWBFE99
5	3 PLY SURGICAL FACEMASK (TIE ON) COLOR: BLUE	IMTBBFE99
6	3 PLY SURGICAL FACEMASK (TIE ON) COLOR: WHITE	IMTWBFE99