



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	Disposable Medical Face Mask
<b>ESTABLISHMENT NAME</b>	ALLIANCE INTEGRITY TRADING SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GA7472622-98671
<b>BRAND NAME</b>	MILKON
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	Disposable medical face mask this product is composed of non-woven cloth, filter cloth, nose clip, mask belt (rubber band); The mask is made of two layers of non-woven cloth and one layer of filter cloth by folding ultrasonic welding. It is non-sterile and disposable. The filtration efficiency of BFE bacteria was $\geq 98\%$ . The filtration efficiency of PFE particles was $\geq 98\%$ .
<b>DEVICE INTENDED PURPOSE</b>	This mask is a disposable product, suitable for the health care of the wearer in the general medical environment without the risk of body fluids and spillage and general health care in public health places.
<b>VALIDITY DATE OF REGISTRATION</b>	27/07/2022 - 26/07/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Disposable Medical Face Mask	MLK-06
2	Disposable Medical Face Mask	MLK-06