

DEVICE DETAILS

NAME OF DEVICE	4 Ply Surgical Face Mask
ESTABLISHMENT NAME	DURIO PPE SDN. BHD.
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA7375622-101645
BRAND NAME	FES
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Surgical mask is fluid-resistant, disposable, and loose-fitting devices that create a physical barrier between the mouth and nose of the wearer and the immediate environment.
DEVICE INTENDED PURPOSE	Disposable device that covers the user's nose and mouth and provides a physical barrier to droplets, fluids, bacteria and particulate materials. Additional non-woven material enhances filtration efficiency & comfort.
VALIDITY DATE OF REGISTRATION	22/08/2022 - 21/08/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	4 Ply Surgical Face Mask - Earloop	4 Ply Surgical Face Mask - Earloop
2	4 Ply Surgical Face Mask - Headloop	4 Ply Surgical Face Mask - Headloop