

DEVICE DETAILS		
NAME OF DEVICE	KF94 Medical Face Mask	
ESTABLISHMENT NAME	ZERO DEGREE HEALTHCARE SDN BHD	
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE	
REGISTRATION NO	GA5573722-90517	
BRAND NAME	MediShield	
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES	
DEVICE GROUPING TYPE	FAMILY	
DEVICE DESCRIPTION	KF94 Medical Face Mask has 4 protective layers with a 3D fit Airy-technology to fit different face shapes. The outer layer protects from the virus while the innermost layer is made of soft material to prevent chafing.	
DEVICE INTENDED PURPOSE	KF94 Medical Face Mask is intended for clinical staff to wear during non-invasive operation. It covers the mouth, nose and mandible of the user, providing a certain physical barrier to prevent the direct penetration of pathogens, microorganisms, particles, etc.	
VALIDITY DATE OF REGISTRATION 18/04/2022 - 17/04/2027		

LIST OF DEVICE		
NO	NAME OF DEVICE DEVICE AUTI	IDENTIFIER
1	KF94 Medical Face Mask	MKF94DG-001
2	KF94 Medical Face Mask	MKF94DB-002
3	KF94 Medical Face Mask	MKF94WR-003
4	KF94 Medical Face Mask	MKF94B-004
5	KF94 Medical Face Mask	MKF94OG-005
6	KF94 Medical Face Mask	MKF94K-006

MKF94G-007

MKF94W-008

MHKF94DG-009

MHKF94DB-010

MHKF94WR-011

Medical Device

KF94 Medical Face Mask

7

8

9

10

11

IST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
12	KF94 Medical Face Mask	MHKF94B-012
13	KF94 Medical Face Mask	MHKF94OG-013
14	KF94 Medical Face Mask	MHKF94K-014
15	KF94 Medical Face Mask	MHKF94W-015

