



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	KF94 4-LAYER PROTECTIVE FACE MASK
<b>ESTABLISHMENT NAME</b>	WATSON'S PERSONAL CARE STORES SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GA6209922-102949
<b>BRAND NAME</b>	WATSONS
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	A single use, non sterile, 4-ply, 3-D face mask with anti-fog design. Made of polypropylene spun bond non-woven fabric and polypropylene meltblown non-woven fabric. Each face mask measuring 185 mm x 102 mm and they are supplied in a pack of 20 pieces.
<b>DEVICE INTENDED PURPOSE</b>	Protecting against bacteria, virus, dust, and micro particles causing allergy.
<b>VALIDITY DATE OF REGISTRATION</b>	06/09/2022 - 05/09/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	WATSONS KF94 4-LAYER PROTECTIVE FACE MASK (NV BLUE)	15163
2	WATSONS KF94 4-LAYER PROTECTIVE FACE MASK (CL BLACK)	15097