



DEVICE DETAILS

NAME OF DEVICE	x'traCool Cooling Gel Patch
ESTABLISHMENT NAME	DURIO PPE SDN. BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA6890722-103958
BRAND NAME	Durio
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	Disposable hydrogel cooling pad that provides cooling effect for up to 8 hours, hydrophilic polymer gel absorbs the energy from the applied fever area, transfer the energy to water molecules and relieves fever via evaporation.
DEVICE INTENDED PURPOSE	Cooling Gel Patch for Relief of Fevers & Headaches
VALIDITY DATE OF REGISTRATION	19/09/2022 - 18/09/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	x'traCool Cooling Gel Patch 	x'traCool