



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	x'traCool Cooling Gel Patch
<b>ESTABLISHMENT NAME</b>	DURIO PPE SDN. BHD.
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GA6890722-103958
<b>BRAND NAME</b>	Durio
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	Disposable hydrogel cooling pad that provides cooling effect for up to 8 hours, hydrophilic polymer gel absorbs the energy from the applied fever area, transfer the energy to water molecules and relieves fever via evaporation.
<b>DEVICE INTENDED PURPOSE</b>	Cooling Gel Patch for Relief of Fevers & Headaches
<b>VALIDITY DATE OF REGISTRATION</b>	19/09/2022 - 18/09/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	x'traCool Cooling Gel Patch	x'traCool