




DEVICE DETAILS

NAME OF DEVICE	Surgical Disposable Face Mask
ESTABLISHMENT NAME	ECOLAND INDUSTRIES SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA5463322-104508
BRAND NAME	mizu
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Surgical Disposable Face Mask is 3-ply protection – (PP waterproff non woven layer/ hight filtration meltblown material/ PP soft absorbent non woven layer); size : 175mm x 95mm
DEVICE INTENDED PURPOSE	Surgical Disposable Face Mask is intended to be worn by operating room personal during surgical procedures as well as dental, isolation and other medical procedures to protect both the patient and the healthcare personnel from transfer of microorganisms, body fluids, and particulate material.
VALIDITY DATE OF REGISTRATION	26/09/2022 - 25/09/2027

LIST OF DEVICE

NO	NAME OF DEVICE 	IDENTIFIER
1	Surgical Disposable Face Mask	MDFM-50-40G
2	Surgical Disposable Face Mask	MDFM-50-40B
3	Surgical Disposable Face Mask	MDFM-50-40W
4	Surgical Disposable Face Mask	MDFMT-50-40G
5	Surgical Disposable Face Mask	MDFMT-50-40B
6	Surgical Disposable Face Mask	MDFMT-50-40W