



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	Surgical Disposable Face Mask
<b>ESTABLISHMENT NAME</b>	ECOLAND INDUSTRIES SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA5463322-104508
<b>BRAND NAME</b>	mizu
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	Surgical Disposable Face Mask is 3-ply protection – ( PP waterproof non woven layer/ high filtration meltblown material/ PP soft absorbent non woven layer); size : 175mm x 95mm
<b>DEVICE INTENDED PURPOSE</b>	Surgical Disposable Face Mask is intended to be worn by operating room personal during surgical procedures as well as dental, isolation and other medical procedures to protect both the patient and the healthcare personnel from transfer of microorganisms, body fluids, and particulate material.
<b>VALIDITY DATE OF REGISTRATION</b>	26/09/2022 - 25/09/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Surgical Disposable Face Mask	MDFM-50-40G
2	Surgical Disposable Face Mask	MDFM-50-40B
3	Surgical Disposable Face Mask	MDFM-50-40W
4	Surgical Disposable Face Mask	MDFMT-50-40G
5	Surgical Disposable Face Mask	MDFMT-50-40B
6	Surgical Disposable Face Mask	MDFMT-50-40W