



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	Silicone Adhesive Tape
<b>ESTABLISHMENT NAME</b>	K-ONE MEDITECH SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GA7228122-104568
<b>BRAND NAME</b>	Huawei
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	It is made of PU film or nonwoven coated with silicone gel. It is applicable to variety of in-dwelling catheters, electrodes or other external fixation that required gentle but robust skin protection and adhesion. Protect the skin which have the risk of ulcer of catheters, or the skin need to be radiation therapy. The tape can be torn into the size you need, as the is cutting line on the film (or tearable nonwoven)
<b>DEVICE INTENDED PURPOSE</b>	Put your scissors away, because with Houseables' soft silicone tape, you won't need them! Every two inches, there is a perforated line so you can easily tear off as much or as little as you may need for the situation. It's quick, easy, and is also shower-ready. Feel free to you through your regular routine with this tape on as it adapts to your skin as you take care of your to-do list.
<b>VALIDITY DATE OF REGISTRATION</b>	26/09/2022 - 25/09/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Silicone Tape Deluxe Type (Skin Color Silicone Tape)	11.05.01
2	Silicone Tape Deluxe Type (Skin Color Silicone Tape)	11.05.02
3	Silicone Tape Deluxe Type (Skin Color Silicone Tape)	11.05.03
4	Silicone Tape Deluxe Type (Skin Color Silicone Tape)	11.05.04
5	Silicone Tape Deluxe Type (Skin Color Silicone Tape)	11.05.05
6	Silicone Tape Deluxe Type (Skin Color Silicone Tape)	11.05.06
7	Silicone Tape Waterproof - Type (PU Film Silicone Tape)	11.06.01
8	Silicone Tape Waterproof - Type (PU Film Silicone Tape)	11.06.02
9	Silicone Tape Waterproof - Type (PU Film Silicone Tape)	11.06.03

<b>NO</b>	<b>NAME OF DEVICE</b>	<b>IDENTIFIER</b>
10	Silicone Tape Waterproof - Type (PU Film Silicone Tape)	11.06.04
11	Silicone Tape Waterproof - Type (PU Film Silicone Tape)	11.06.05
12	Silicone Tape Waterproof - Type (PU Film Silicone Tape)	11.06.06
13	Silicone Tape Economic - Type (Non Woven Silicone Tape)	11.07.01
14	Silicone Tape Economic - Type (Non Woven Silicone Tape)	11.07.02
15	Silicone Tape Economic - Type (Non Woven Silicone Tape)	11.07.03
16	Silicone Tape Economic - Type (Non Woven Silicone Tape)	11.07.04
17	Silicone Tape Economic - Type (Non Woven Silicone Tape)	11.07.05
18	Silicone Tape Economic - Type (Non Woven Silicone Tape)	11.07.06



**MEDICAL DEVICE AUTHORITY**