



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	4-PLY PREMIUM MEDICAL DISPOSABLE FACE MASK
<b>ESTABLISHMENT NAME</b>	JL WATER ENGINEERING SDN. BHD.
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA5340522-107407
<b>BRAND NAME</b>	GENOVIQ
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	A four-layered loose fitting mask that is designed to cover wearer's nose and mouth that meets fluid barrier or filtration efficiency levels
<b>DEVICE INTENDED PURPOSE</b>	Used as barrier for user working in general medical environment to avoid unwanted inhalation or protecting from droplet and spill to avoid any unexpected infection of flu or disease. This product is disposable and suitable for health care of the wearer in the general medical environment and the general care in public health places where there is any risk of bodily fluids and droplets.
<b>VALIDITY DATE OF REGISTRATION</b>	08/11/2022 - 07/11/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	4-ply Premium Medical Disposable Face Mask	BBK01
2	4-ply Premium Medical Disposable Face Mask	BWH02
3	4-ply Premium Medical Disposable Face Mask	BPK03
4	4-ply Premium Medical Disposable Face Mask	BCY04
5	4-ply Premium Medical Disposable Face Mask	BPL05
6	4-ply Premium Medical Disposable Face Mask	CBO01
7	4-ply Premium Medical Disposable Face Mask	CWP02
8	4-ply Premium Medical Disposable Face Mask	CPB03
9	4-ply Premium Medical Disposable Face Mask	CTG04
10	4-ply Premium Medical Disposable Face Mask	CPG05