




DEVICE DETAILS

NAME OF DEVICE	Proxima 3Ply Disposable Medical/Surgical Face Mask
ESTABLISHMENT NAME	PROXIMA GLOBAL SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA7465922-107979
BRAND NAME	Proxima
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	The disposal medical/surgical face mask consist of non-woven fabric and electrostatic melt blown filter, along with elastic ear loop and adjustable nose piece to fit the nose.
DEVICE INTENDED PURPOSE	The intended use of medical/ surgical face mask to protect general healthcare, patient and public users against of airborne particulates, microorganisms, blood and body fluid to reduce the risk of infections.
VALIDITY DATE OF REGISTRATION	08/11/2022 - 07/11/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Proxima 3Ply Disposable Medical/Surgical Face Mask 	Proxima 3Ply Surgical face mask