



DEVICE DETAILS

NAME OF DEVICE	FINGERTIP PULSE OXIMETER
ESTABLISHMENT NAME	KINSMEDIC SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GB2362522-109761
BRAND NAME	CHOICEMMED
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	A pulse oximeter uses two frequencies of light (red and infrared) to determine the percentage (%) of hemoglobin in the blood that is saturated with oxygen. The percentage is called blood oxygen saturation, or SpO ₂ . A pulse oximeter also measures and displays the pulse rate at the same time it measures the SpO ₂ level.
DEVICE INTENDED PURPOSE	Fingertip pulse oximeter is a handheld non-invasive device intended for spot checking of oxygen saturation of arterial hemoglobin (SpO ₂) and pulse rate of adult and pediatric patients in hospitals, hospital type facilities and homecare.
VALIDITY DATE OF REGISTRATION	22/11/2022 - 21/11/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Fingertip Pulse Oximeter	MD300CN310
2	Fingertip Pulse Oximeter	MD300C52
3	Children Fingertip Pulse Oximeter	MD300C55
4	Pulse Oximeter	MD300C15D
5	Pulse Oximeter	MD300C29