



DEVICE DETAILS

| | |
|--------------------------------------|---|
| NAME OF DEVICE | FINGERTIP PULSE OXIMETER |
| ESTABLISHMENT NAME | KINSMEDIC SDN BHD |
| ROLE OF ESTABLISHMENT | AUTHORIZED REPRESENTATIVE |
| REGISTRATION NO | GB2362522-109761 |
| BRAND NAME | CHOICEMMED |
| MEDICAL DEVICE CATEGORY | MD 1100 - GENERAL ACTIVE MEDICAL DEVICES |
| DEVICE GROUPING TYPE | FAMILY |
| DEVICE DESCRIPTION | A pulse oximeter uses two frequencies of light (red and infrared) to determine the percentage (%) of hemoglobin in the blood that is saturated with oxygen. The percentage is called blood oxygen saturation, or SpO ₂ . A pulse oximeter also measures and displays the pulse rate at the same time it measures the SpO ₂ level. |
| DEVICE INTENDED PURPOSE | Fingertip pulse oximeter is a handheld non-invasive device intended for spot checking of oxygen saturation of arterial hemoglobin (SpO ₂) and pulse rate of adult and pediatric patients in hospitals, hospital type facilities and homecare. |
| VALIDITY DATE OF REGISTRATION | 22/11/2022 - 21/11/2027 |

LIST OF DEVICE

| NO | NAME OF DEVICE | IDENTIFIER |
|----|-----------------------------------|------------|
| 1 | Fingertip Pulse Oximeter | MD300CN310 |
| 2 | Fingertip Pulse Oximeter | MD300C52 |
| 3 | Children Fingertip Pulse Oximeter | MD300C55 |
| 4 | Pulse Oximeter | MD300C15D |
| 5 | Pulse Oximeter | MD300C29 |