

DEVICE DETAILS		
NAME OF DEVICE	FINGERTIP PULSE OXIMETER	
ESTABLISHMENT NAME	KINSMEDIC SDN BHD	
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE	
REGISTRATION NO	GB2362522-109761	
BRAND NAME	CHOICEMMED	
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES	
DEVICE GROUPING TYPE	FAMILY	
DEVICE DESCRIPTION	A pulse oximeter uses two frequencies of light (red and infrared) to determine the percentage (%) of hemoglobin in the blood that is saturated with oxygen. The percentage is called blood oxygen saturation, or SpO2. A pulse oximeter also measures and displays the pulse rate at the same time it measures the SpO2 level.	
DEVICE INTENDED PURPOSE	Fingertip pulse oximeter is a handheld non-invasive device intended for spot checking of oxygen saturation of arterial hemoglobin (SpO2) and pulse rate of adult and pediatric patients in hospitals, hospital type facilities and homecare.	
VALIDITY DATE OF REGISTRATION	22/11/2022 - 21/11/2027	

LIST OF DEVICE		
NO	NAME OF DEVICE	IDENTIFIER
1	Fingertip Pulse Oximeter	MD300CN310
2	Fingertip Pulse Oximeter	MD300C52
3	Children Fingertip Pulse Oximeter	MD300C55
4	Pulse Oximeter	MD300C15D
5	Pulse Oximeter	MD300C29