



## DEVICE DETAILS

NAME OF DEVICE	4 Layer Surgical Face Mask
ESTABLISHMENT NAME	PROXIMA GLOBAL SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA3496022-109952
BRAND NAME	Proxima
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	The disposal medical/surgical face mask consist of non-woven fabric and electrostatic melt blown filter, along with elastic ear loop and adjustable nose piece to fit the nose.
DEVICE INTENDED PURPOSE	The intended use of medical/ surgical face mask to protect general healthcare, patient and public users against of airborne particulates, microorganisms, blood and body fluid to reduce the risk of infections.
VALIDITY DATE OF REGISTRATION	28/11/2022 - 27/11/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	4 Layer Surgical Face Mask	 4ply Disposable Surgical face mask