



DEVICE DETAILS

NAME OF DEVICE	KF94 PREMIUM FACE MASK
ESTABLISHMENT NAME	BEWONDER HOLDING SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA2416022-110606
BRAND NAME	BEMEDIC
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	KF94 Face Mask. The outer layer with high quality non-woven fabric that repels fluid for a higher level of protection. Second and third middle layers with melt-blown fabric for effective particle filtration. Inner layer with extra soft premium non-woven fabric that absorbs moisture and helps reduce skin irritation from long wear.
DEVICE INTENDED PURPOSE	Use for bacterial filtration with submicron protection against airborne particles, viruses and bacteria. Designed to keep out biofluid and other droplets.
VALIDITY DATE OF REGISTRATION	07/12/2022 - 06/12/2027

LIST OF DEVICE



NO	NAME OF DEVICE	IDENTIFIER
1	BEMEDIC KF94 KOREAN PREMIUM FACE MASK	12122119
2	BEMEDIC KF94 HIJAB HEADLOOP KOREAN PREMIUM FACE MASK	12122120
3	BEMEDIC KF94 DISPOSABLE FACE MASK	12121949
4	BEMEDIC KF94 DISPOSABLE FACE MASK	12121947
5	BEMEDIC KF94 DISPOSABLE FACE MASK	12121948
6	BEMEDIC KF94 DISPOSABLE FACE MASK	12121905
7	BEMEDIC KF94 DISPOSABLE FACE MASK	12121907
8	BEMEDIC KF94 DISPOSABLE FACE MASK	12121946
9	BEMEDIC KF94 DISPOSABLE FACE MASK	12121906
10	BEMEDIC KF94 DISPOSABLE FACE MASK	12121945
11	BEMEDIC KF94 DISPOSABLE FACE MASK	12121945
12	BEMEDIC KF94 KOREAN PREMIUM FACE MASK - LIMITED EDITION	12122125

NO	NAME OF DEVICE	IDENTIFIER
13	BEMEDIC KF94 KOREAN PREMIUM FACE MASK - LIMITED EDITION	12122123
14	BEMEDIC KF94 KOREAN PREMIUM FACE MASK - LIMITED EDITION	12122127
15	BEMEDIC KF94 KOREAN PREMIUM FACE MASK - LIMITED EDITION	12122126
16	BEMEDIC KF94 KOREAN PREMIUM FACE MASK - LIMITED EDITION	12122124
17	BEMEDIC KF94 KOREAN PREMIUM FACE MASK - LIMITED EDITION	12122130
18	BEMEDIC KF94 KOREAN PREMIUM FACE MASK - LIMITED EDITION	12122128
19	BEMEDIC KF94 KOREAN PREMIUM FACE MASK - LIMITED EDITION	12122132
20	BEMEDIC KF94 KOREAN PREMIUM FACE MASK - LIMITED EDITION	12122131
21	BEMEDIC KF94 KOREAN PREMIUM FACE MASK - LIMITED EDITION	12122129
22	BEMEDIC LUNAR NEW YEAR KF94 - LIMITED EDITION	12122169
23	BEMEDIC LUNAR NEW YEAR KF94 - LIMITED EDITION	12122168
24	BEMEDIC LUNAR NEW YEAR KF94 - LIMITED EDITION	12122170
25	BEMEDIC LUNAR NEW YEAR KF94 - LIMITED EDITION	12122167
26	BEMEDIC LUNAR NEW YEAR KF94 - LIMITED EDITION	12122173
27	BEMEDIC LUNAR NEW YEAR KF94 - LIMITED EDITION	12122172
28	BEMEDIC LUNAR NEW YEAR KF94 - LIMITED EDITION	12122174
29	BEMEDIC LUNAR NEW YEAR KF94 - LIMITED EDITION	12122171
30	BEMEDIC KF94 KOREAN PREMIUM FACE MASK	12122155
31	BEMEDIC KF94 KOREAN PREMIUM FACE MASK	12122184