




DEVICE DETAILS

NAME OF DEVICE	Adult 3-Ply Premium Medical Face Mask
ESTABLISHMENT NAME	H&A PRODUCTS SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA9574122-111543
BRAND NAME	H&A Products
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	H&A Medical Face Mask is made of 3 layers of high quality non-woven fabric , meeting the requirements of ASTM F2100 and BS EN 14683. It achieves Bacterial Filtration Efficiency (BFE) of more than 99% and Sub-micron Particulate Filtration Efficiency (PFE) of more than 98%
DEVICE INTENDED PURPOSE	H&A Adult 3-Ply Premium Medical Face Mask is designed for single use and it provides optimum protection against transfer of microorganisms, blood fluid and airborne particulate.
VALIDITY DATE OF REGISTRATION	14/12/2022 - 13/12/2027

LIST OF DEVICE

NO	NAME OF DEVICE 	IDENTIFIER
1	Adult 3-Ply Premium Medical Face Mask	FMFF3L02AH