



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	Adult 3-Ply Premium Medical Face Mask
<b>ESTABLISHMENT NAME</b>	H&A PRODUCTS SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA9574122-111543
<b>BRAND NAME</b>	H&A Products
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	H&A Medical Face Mask is made of 3 layers of high quality non-woven fabric , meeting the requirements of ASTM F2100 and BS EN 14683. It achieves Bacterial Filtration Efficiency (BFE) of more than 99% and Sub-micron Particulate Filtration Efficiency (PFE) of more than 98%
<b>DEVICE INTENDED PURPOSE</b>	H&A Adult 3-Ply Premium Medical Face Mask is designed for single use and it provides optimum protection against transfer of microorganisms, blood fluid and airborne particulate.
<b>VALIDITY DATE OF REGISTRATION</b>	14/12/2022 - 13/12/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Adult 3-Ply Premium Medical Face Mask	FMFF3L02AH