



# Public Search



Malaysia Medical Device Register (MMDR)

## DEVICE DETAILS

NAME OF DEVICE	KF94 Face Mask
ESTABLISHMENT NAME	TST HEALTH CARE SDN. BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA8868622-111713
BRAND NAME	DYMEX
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	Medical Face Mask - KF94 are disposable medical barrier made up of 4-ply fluid resistance outer layer, excellent filtering of non woven and melt blown fabric that provides effective protection and elastic ear-loop for extra comfort of the users. This mask is intended to protect wearer's mouth and nose from any infectious air-borne viruses, bacteria, dust or any other air-borne particles.
DEVICE INTENDED PURPOSE	This mask is intended to protect wearer's mouth and nose from any infectious air-borne viruses, bacteria, dust or any other air-borne particles.
VALIDITY DATE OF REGISTRATION	14/12/2022 - 13/12/2027

## LIST OF DEVICE



NO	NAME OF DEVICE	IDENTIFIER
1	KF94 Face Mask	MLKF94-W