



DEVICE DETAILS

NAME OF DEVICE	HIMAYA 3-PLY PREMIUM HIJAB MEDICAL FACE MASK
ESTABLISHMENT NAME	NEUTROVIS SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA7484122-113389
BRAND NAME	HIMAYA
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Himaya 3ply Premium Hijab Medical Face Mask is made of premium medical grade melt-blown filter fabric for better breathability, comfort and reduced irritation
DEVICE INTENDED PURPOSE	Himaya 3ply Premium Hijab Medical Face Mask with BFE $\geq 98\%$ and PFE \geq is able to filter out live bacteria (viable particles) from 1 to 5 microns and micro particles (nonviable particles) from 0.1 micron to 1 micron
VALIDITY DATE OF REGISTRATION	27/12/2022 - 26/12/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Himaya 3-Ply Premium Hijab Medical Face Mask	HMY001
2	Himaya 3-Ply Premium Hijab Medical Face Mask	HMY002
3	Himaya 3-Ply Premium Hijab Medical Face Mask	HMY003
4	Himaya 3-Ply Premium Hijab Medical Face Mask	HMY004
5	Himaya 3-Ply Premium Hijab Medical Face Mask	HMY005
6	Himaya 3-Ply Premium Hijab Medical Face Mask	HMY006
7	Himaya 3-Ply Premium Hijab Medical Face Mask	HMY007
8	Himaya 3-Ply Premium Hijab Medical Face Mask	HMY008
9	Himaya 3-Ply Premium Hijab Medical Face Mask	HMY009
10	Himaya 3-Ply Premium Hijab Medical Face Mask	HMY010
11	Himaya 3-Ply Premium Hijab Medical Face Mask	HMY014