



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	NEUTROVIS 3-PLY PREMIUM HIJAB MEDICAL FACE MASK
<b>ESTABLISHMENT NAME</b>	NEUTROVIS SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA7835622-113407
<b>BRAND NAME</b>	NEUTROVIS
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	Neutrovis 3ply Premium Hijab Medical Face Mask is made of premium medical grade melt-blown filter fabric for better breathability, comfort and reduced irritation
<b>DEVICE INTENDED PURPOSE</b>	Neutrovis 3ply Premium Hijab Medical Face Mask with BFE $\geq$ 98% and PFE $\geq$ 96% is able to filter out live bacteria (viable particles) from 1 to 5 microns and micro particles (nonviable particles) from 0.1 micron to 1 micron
<b>VALIDITY DATE OF REGISTRATION</b>	27/12/2022 - 26/12/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Neutrovis 3-Ply Premium Hijab Medical Face Mask	NTV16
2	Neutrovis 3-Ply Premium Hijab Medical Face Mask	NTV24
3	Neutrovis 3-Ply Premium Hijab Medical Face Mask	NTV94