



DEVICE DETAILS

NAME OF DEVICE	NEUTROVIS 3-PLY PREMIUM HIJAB MEDICAL FACE MASK
ESTABLISHMENT NAME	NEUTROVIS SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA7835622-113407
BRAND NAME	NEUTROVIS
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Neutrovis 3ply Premium Hijab Medical Face Mask is made of premium medical grade melt-blown filter fabric for better breathability, comfort and reduced irritation
DEVICE INTENDED PURPOSE	Neutrovis 3ply Premium Hijab Medical Face Mask with BFE \geq 98% and PFE \geq 96% is able to filter out live bacteria (viable particles) from 1 to 5 microns and micro particles (nonviable particles) from 0.1 micron to 1 micron
VALIDITY DATE OF REGISTRATION	27/12/2022 - 26/12/2027

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Neutrovis 3-Ply Premium Hijab Medical Face Mask	NTV16
2	Neutrovis 3-Ply Premium Hijab Medical Face Mask	NTV24
3	Neutrovis 3-Ply Premium Hijab Medical Face Mask	NTV94